Defibrillator Analyzer

DELTA 3000



DELTA 3000, Defibrillator Analyzer / Pacemaker Tester is the latest addition of our Defibrillator Analyzer with full Pacer Tester and ECG Simulator built in.

DELTA 3000, Defibrillator Analyzer / Pacemaker Tester combines the complete testing functions of a Defibrillator Tester, a Pacemaker Analyzer, and a 12 Lead ECG / Arrhythmia simulator into one compact, lightweight, and easy to use instrument. No other instrument delivers the performance testing power of the **DELTA 3000** Defibrillator Analyzer / Pacemaker Tester.

DELTA 3000, Defibrillator Analyzer / Pacemaker Tester tests all defibrillators including AEDs. It measures both monophasic and bi-phasic waveforms. It tests all external pacemakers including Transthoracic, Transvenous, and AV sequential models. It has 12 built-in user selectable test loads for Transthoracic pacing testing. The built-in ECG / Arrhythmia simulator allows performance testing of any ECG equipment.

Features:

Compact, Lightweight, and Rugged Simple, Menu Driven Operation Large 4-line LCD Simultaneous Display of All Test Measurements Dedicated Print Key for Immediate Documentation RS 232 Interface Battery or AC Operation Two-Year Warranty



Measurements:

Measures: Energy, Peak Voltage, Peak Current, and Cardio Sync. Time Storage and Playback of output waveforms Tests all External Pacemakers - Transthorasic, Transvenous, and AV Pacers Measures: Amplitude, Rate, Pulse Width, Pulse Energy, Sensed and Paced Refractory Periods 50 and 60 Hz Immunity Tests Transthorasic Test Loads: Fixed 50 Ohm or 12 User Selectable Loads to 1000 Ohms A / V Pacer Test Load: 500 Ohm ECG Simulator: Tests all ECG Equipment, 12 lead ECG, User Selectable Rates and Amplitudes Performance Waveforms: Sine, Square, Triangle, and Pulse, 12 Arrhythmias



深圳市展业达鸿科技有限公司

吕先生:15920060912(微信同号) 0755-22934005(座机) 地址:深圳市福田区八卦二路八卦岭工业区615栋419 邮箱:hongqi@thingstet.com 网址:www.thingstest.com

Defibrillator Analyzer

DELTA 3000

Specifications:

Defibrillator Analyzer

Energy Measurement: Load: 50 Ohm ± 1 % Range (High): 900 J, Res.1 J Resolution: 1 Joules, Range (Low): 0 - 100 Resolution: 0.1 Joules, Accuracy: $\pm 2\%$ of the reading + 1 J for Hi Range or + 0.5 J for Low Range Maximum Voltage: 5200 V Maximum Current: 100 Amp (H) Sync. Time: 0 - 250 ms Waveform Output: Real Time: Output 1000 : 1 Stored. Discharge waveform expanded 200:1 and is repeated 8 sec. It can be viewed through the paddles and through the ECG High Level Output

Pacemaker Analyzer

Transcutaneous: Test Load: 50 Ohm (Fixed) and variable 100 to 1000 Ohms user selectable Current: 1 to 200 mA, Accuracy: $\pm 5\% + 0.5$ mA Pulse Width: 0.5 to 80 ms, Accuracy: $\pm 1\% + 0.1$ ms Refractory: 20 - 500 ms, Accuracy: $\pm 2\%$ or ± 10 ms Atrial and Ventricular: Test Load: 500 Ohms Current: 1 to 25 mA, Accuracy: \pm 5% + 0.5 mA Pulse Rate: 30 to 800 ppm, Accuracy: $\pm 1\% + 0.1$ ms and Pulse Refractory: ± 20 - 500 ms, Accuracy: $\pm 2\%$ or 10 ms Immunity Test: 50 / 60 Hz Arrhythmias:

General

Power Requirements: Two-9V Alkaline Batteries Or AC Adapter Dimension: 10 x 11 x 2.25 inches (25.4 x 27.9 x 5.7 cm) Weight: 4 lbs (1.8 kg

ECGArrhythmia Simulator

ECG: Normal Sinus Rhythm. Rate: 30, 60, 90, 120, 240 Performance: Sine, Square, Triangle Rate: 0.5, 1, 2, 10, 15, 20, 25, 40 Hz VFBC, VFBF, VTAC, AFIB, PVC, PVC1, PVC2, BGY, RBB, AFIB, ATRFT, NSR Waveform Output: 12 Lead ECG and Paddle and High Level Output

Environmental:

Operating Temperature: 59° to 95° F (15° to 35° C), Storage Temperature: 32° to 131° F (0° to 55° C)

Ordering Information:

Description

Defibrillator / Pacemaker Analyzer with Transcutaneous A-V Pacer Analyzer Defibrillator / Pacemaker Analyzer with Transcutaneous/Transvenous Analyzer



深圳市展业达鸿科技有限公司

吕先生:15920060912(微信同号) 0755-22934005(座机) 地址:深圳市福田区八卦二路八卦岭工业区615栋419 邮箱:hongqi@thingstet.com 网址:www.thingstest.com DELAIA 3000 Page2 of 2